

USP-100SP Features:

- Digital reading
- Continuous or Pulsed Power
- · Integrated carrying case
- 0.2 to 30.0 Watt Range
- Accuracy = ± 3% + one count reading
- 0.5 to 10 MHz Frequency Response
- 3" Max Transducer Diameter
- Built-in transducer holder
- Degassed water is the testing medium
- Line Power or Battery Operation
- Compact size of 8.5" H x 12" L x 8.75" W
- Unit weight is 11 lbs

USP-30 Features:

- Analog reading
- 0.15 to 30.0 Watt Range
- Continuous or Pulsed Power
- 1 to 10 MHz Frequency Response
- Integrated carrying case
- Built-in transducer holder

ตัวแทนจำหน่าย

Perfect Quality (Thailand) Co., Ltd. 99 Moo. 11 Dongkeng Subdistrict, Nong Song Hong District, Khon Kaen Provine, 40190 Phone: 043-041451 Email: Perfectquality_@hotmail.com www. สอบเทียบเครื่องมือแพทย์.com

Ultrasound Power Meters

BC Biomedical USP-100SP & USP-30

The USP-100SP is designed to measure the ultrasound power output of diagnostic or therapeutic transducers up to 30 watts and according to AIUM, JCAHO, FDA and NIST guide-lines. The USP-100SP has a display resolution of 50mW. The compact and light weight USP-100SP meets all the requirements of a portable test instrument with reliability and proven performance.

The USP-30 is designed to measure the ultrasonic power output of diagnostic or therapeutic transducers up to 30 Watts using the radiation force balance method. The USP- 30 consists of a precision mechanical balance, which measures from 0.01 to 10.00 grams, a conical target suspended on a triangular hanger, a water tank, a transducer support rod w/clamp and a carrying case. The transducer to be tested is center mounted over the conical target with the transducer's radiating area immersed facedown in the water medium. This system provides a sturdy instrument with good measurement accuracy at an economical cost.





Specifications

USP-100SP

Measure Range	0.2 to 30 Watts
Resolution	± 50 mW
Accuracy	± 3 + one count reading
Transducer Frequency	
Maximum Transducer Diameter	3"
Test Medium	Degassed, destilled, room temp water
Stabilization Time	2.5 seconds
Display	3 digit LCD, 1.3" high
Line Power	120 or 240 VAC plug-in adapter
Battery	4 AA, 20 hrs. continuous battery life
Dimensions	8.5" H x 12" L x 8.75" W
Weight	11 lbs.

USP-30

Power Range	150 mW to30 Watts
Resolution	. ±150 mW
Minimum Detectable Power Level	±150 mW
Minimum Accuracy	. ±5%
Expected Accuracy	
Display	. Vernier Dial
Display Sensitivity	
Zeroing Method	Manually adjusted
Test Tank Lining	0.5" neoprene rubber
Test Medium	
	Degassed, destilled, room temp water
Test Medium	Degassed, destilled, room temp water
Test Medium Test Tank Liquid Capacity	Degassed, destilled, room temp water 1350 mL
Test Medium Test Tank Liquid Capacity Test Tank Dimensions	Degassed, destilled, room temp water1350 mL5.75" dia. x 5.55"
Test Medium Test Tank Liquid Capacity Test Tank Dimensions Outside Inside	Degassed, destilled, room temp water1350 mL5.75" dia. x 5.55"4.5" dia x 5.25 D.
Test Medium Test Tank Liquid Capacity Test Tank Dimensions Outside Inside Cone Target Dimensions	Degassed, destilled, room temp water 1350 mL 5.75" dia. x 5.55" 4.5" dia x 5.25 D. 3.5" dia. x 2" H
Test Medium Test Tank Liquid Capacity Test Tank Dimensions Outside Inside	Degassed, destilled, room temp water 1350 mL 5.75" dia. x 5.55" 4.5" dia x 5.25 D. 3.5" dia. x 2" H 19" H x 17" L x 7" W

About BC Group International

BC Group is a leader in manufacturing, servicing and selling biomedical testing equipment. With dealers across the world, BC Group has achieved international success and continues to grow. Our equipment is designed to thoroughly test medical equipment and maximize your investments. We have a promise to provide quality at all levels of operation and help ensure equipment is safe for patients and medical staff.

BC Group is committed to exceeding customer expectations and enhancing customer satisfaction through the continual improvement of our products, services and Quality Management System.

ตัวแทนจำหน่าย

Perfect Quality (Thailand) Co., Ltd. 99 Moo. 11 Dongkeng Subdistrict, Nong Song Hong District, Khon Kaen Provine, 40190 Phone: 043-041451 Email: Perfectquality_@hotmail.com www. สอบเทียบเครื่องมือแพทย์.com